

PATCH INFORMATION

- The majority of the infusion happens in the first 5-15 minutes on each treatment area, but the patches may remain on the face for up to 60 minutes.
- Retinol easily breaks down. Please store your patches in a cool dry place away from sunlight. You may even prefer to store them in the refrigerator to have a cooling gel patch for treatment.
- Retinol can make your skin more sensitive to the sun. After use, we recommend using a high factor sunscreen.
- Ensure patches are well adhered to the face before treatment so they don't peel off.
- The patches are designed for single use only.

GENERAL INFORMATION

- Clean by wiping with a damp cloth. Do NOT use solvents or detergents.
- Do NOT immerse the unit in water.
- Do NOT dissemble or repair. There are no user serviceable parts inside the unit. If you have any problems using the product, contact our Customer Care Department.
- If the unit is damaged do not use it.
- Do NOT use near to shortwave or microwave equipment.
- If you experience unexpected behaviour of this product then it is possible some external electromagnetic interference is responsible e.g. electric motors/ loudspeakers. You should remove this source and restart the product to restore full operation.
- For best results we recommend completing the cleansing pre-treatment before moving onto the infusion treatment.

BATTERY INFORMATION



To fit the battery, release the clip to open the battery cover. Insert 2 AAA alkaline batteries ensuring the (+) terminal is as shown. Replace cover.

- If the unit is left unused for a long period of time, remove the batteries.
- Do NOT use rechargeable batteries.
- Please ensure that the batteries are disposed of in the correct manner in accordance with local authority requirements.

PATCH INFORMATION

- The majority of the infusion happens in the first 5-15 minutes on each treatment area, but the patches may remain on the face for up to 60 minutes.
- Retinol easily breaks down. Please store your patches in a cool dry place away from sunlight. You may even prefer to store them in the refrigerator to have a cooling gel patch for treatment.
- Retinol can make your skin more sensitive to the sun. After use, we recommend using a high factor sunscreen.
- Ensure patches are well adhered to the face before treatment so they don't peel off.
- The patches are designed for single use only.

GENERAL INFORMATION

- Clean by wiping with a damp cloth. Do NOT use solvents or detergents.
- Do NOT immerse the unit in water.
- Do NOT dissemble or repair. There are no user serviceable parts inside the unit. If you have any problems using the product, contact our Customer Care Department.
- If the unit is damaged do not use it.
- Do NOT use near to shortwave or microwave equipment.
- If you experience unexpected behaviour of this product then it is possible some external electromagnetic interference is responsible e.g. electric motors/ loudspeakers. You should remove this source and restart the product to restore full operation.
- For best results we recommend completing the cleansing pre-treatment before moving onto the infusion treatment.

BATTERY INFORMATION



To fit the battery, release the clip to open the battery cover. Insert 2 AAA alkaline batteries ensuring the (+) terminal is as shown. Replace cover.

- If the unit is left unused for a long period of time, remove the batteries.
- Do NOT use rechargeable batteries.
- Please ensure that the batteries are disposed of in the correct manner in accordance with local authority requirements.

CUSTOMER CARE

If you have any questions, or if you experience any difficulties please call our Care line on: **+44 (0)1242 702345** or email us at: **online@riobeauty.com**

To register your two year guarantee online visit: **www.riobeauty.com/warranty_online.htm**

DISPOSAL



The Waste of Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill.

The product is classed as Electrical or Electronic equipment so please ensure that at the end of the life of this product it is disposed of in the correct manner in accordance with local authority requirements. It must NOT be disposed of with household waste.



Please refer to these instructions.




Rio, Dezac and Galvanic Wand are trademarks of The Dezac Group Ltd
© The Dezac Group Ltd 2012
PO Box 17, Cheltenham Spa, England GL53 7ET

www.riobeauty.com


1-IN-GVWD-UK/1.1

GALVANIC WAND™



The Galvanic Wand uses gentle micro currents to infuse Collagen and Retinol into the deeper layers of the skin, for an intense wrinkle reduction and skin rejuvenation treatment.

GALVANIC WAND™



The Galvanic Wand uses gentle micro currents to infuse Collagen and Retinol into the deeper layers of the skin, for an intense wrinkle reduction and skin rejuvenation treatment.



The Waste of Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill.

The product is classed as Electrical or Electronic equipment so please ensure that at the end of the life of this product it is disposed of in the correct manner in accordance with local authority requirements. It must NOT be disposed of with household waste.



Please refer to these instructions.



Rio, Dezac and Galvanic Wand are trademarks of The Dezac Group Ltd
© The Dezac Group Ltd 2012
PO Box 17, Cheltenham Spa, England GL53 7ET

www.riobeauty.com

1-IN-GVWD-UK/1.1

INTRODUCTION

The Galvanic Wand treatment regime comprises two distinct treatments.


- Iontophoresis Cleansing
- Iontophoresis Infusion.

The cleansing treatment helps remove impurities from the skin and is also the perfect preparation for the infusion treatment.

Iontophoresis Infusion uses small micro currents to carry the active ingredients for skin repair and rejuvenation deep into the skin. This transdermal infusion greatly increase the effectiveness of topically applied treatments which would otherwise find it difficult to penetrate the skin.

Short-term infusion results can usually be seen almost immediately. With regular use the effects can become long lasting.

TREATMENT ZONES



IMPORTANT - PATCH TEST
Before you start you must ALWAYS do a patch test on an inconspicuous area of your skin for possible skin reaction. Wait 24 hours. If irritation occurs do not continue treatment.

INTRODUCTION

The Galvanic Wand treatment regime comprises two distinct treatments.


- Iontophoresis Cleansing
- Iontophoresis Infusion.

The cleansing treatment helps remove impurities from the skin and is also the perfect preparation for the infusion treatment.

Iontophoresis Infusion uses small micro currents to carry the active ingredients for skin repair and rejuvenation deep into the skin. This transdermal infusion greatly increase the effectiveness of topically applied treatments which would otherwise find it difficult to penetrate the skin.


Short-term infusion results can usually be seen almost immediately. With regular use the effects can become long lasting.

TREATMENT ZONES



IMPORTANT - PATCH TEST
Before you start you must ALWAYS do a patch test on an inconspicuous area of your skin for possible skin reaction. Wait 24 hours. If irritation occurs do not continue treatment.


KNOW YOUR UNIT



PRECAUTIONS

- **Please read all instructions thoroughly and retain for future reference.**
- Do NOT use if you are pregnant or if you have any of the following: pacemaker, heart disorder, receiving treatment for low blood pressure, epilepsy, multiple sclerosis, or if you are recovering from an operation.
- Do NOT use on areas of skin with disorders such as rosacea, pustular acne, abraded skin, sensitive skin, or over scars, moles or varicose veins.
- If the gel or patches come into contact with the eyes, rinse thoroughly with water. If irritation persists, consult a doctor.
- Keep out of reach of children and pets.
- This appliance is not intended for use by persons, including children, with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instructions concerning use of the appliance by a person responsible for their safety.
- Take particular care when moving near the eye. Ensure the roller is at least 1cm away from the eye at all times.

KNOW YOUR UNIT



PRECAUTIONS


- **Please read all instructions thoroughly and retain for future reference.**
- Do NOT use if you are pregnant or if you have any of the following: pacemaker, heart disorder, receiving treatment for low blood pressure, epilepsy, multiple sclerosis, or if you are recovering from an operation.
- Do NOT use on areas of skin with disorders such as rosacea, pustular acne, abraded skin, sensitive skin, or over scars, moles or varicose veins.
- If the gel or patches come into contact with the eyes, rinse thoroughly with water. If irritation persists, consult a doctor.
- Keep out of reach of children and pets.
- This appliance is not intended for use by persons, including children, with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instructions concerning use of the appliance by a person responsible for their safety.
- Take particular care when moving near the eye. Ensure the roller is at least 1cm away from the eye at all times.

CLEANSING PRE-TREATMENT

We recommend the Vitamin Cleansing Gel is used only once in 24 hours. If you have sensitive skin we would recommend leaving 48 hours between treatments.



STEP 1 Onto clean skin apply a generous amount of Vitamin Cleansing Gel to your face.
NOTE: Without the gel the treatment may feel prickly.



STEP 2 Push the power button to ON. Set the micro-current button to (-) negative and adjust the intensity wheel to level 1.



STEP 3 Move the roller slowly back and forth over the treatment area.
NOTE: It is important that you hold the unit by the conductive panels. Correct operation is indicated by a red light when the roller contacts the skin/gel. The unit will not make a noise or vibrate. If insufficient gel is used the treatment may feel prickly.



STEP 4 To increase the intensity turn the intensity wheel clockwise.
Working over the whole face should take approximately 3-6 minutes.




STEP 5 Thoroughly clean and dry skin. You can continue directly onto infusion treatment.

CLEANSING PRE-TREATMENT

We recommend the Vitamin Cleansing Gel is used only once in 24 hours. If you have sensitive skin we would recommend leaving 48 hours between treatments.



STEP 1 Onto clean skin apply a generous amount of Vitamin Cleansing Gel to your face.
NOTE: Without the gel the treatment may feel prickly.



STEP 2 Push the power button to ON. Set the micro-current button to (-) negative and adjust the intensity wheel to level 1.



STEP 3 Move the roller slowly back and forth over the treatment area.
NOTE: It is important that you hold the unit by the conductive panels. Correct operation is indicated by a red light when the roller contacts the skin/gel. The unit will not make a noise or vibrate. If insufficient gel is used the treatment may feel prickly.




STEP 4 To increase the intensity turn the intensity wheel clockwise.
Working over the whole face should take approximately 3-6 minutes.




STEP 5 Thoroughly clean and dry skin. You can continue directly onto infusion treatment.

INFUSION TREATMENT

We recommend Collagen and Retinol patches are only used once in 24 hours. If you have sensitive skin we would recommend leaving 48 hours between treatments.



STEP 1 Tear along the perforated edge of the Collagen and Retinol patch sachet. Take a patch from the plastic casing and peel off the liners to reveal the gel patch.



STEP 2 Place the patch onto the area to be treated. Place the other patch on the opposite side of the face.



STEP 3 Push the power button to ON. Set the micro-current button to (+)positive and set the intensity wheel to level 1



STEP 4 Slowly and gently move the roller back and forth on the patch.

Repeat the process over the other patch. Continue treatment by alternating between patches for 10-30 minutes (5-15 minutes on each patch).

NOTE: It is important that you hold the unit by the conductive panels. Correct operation is indicated by a blue light when the roller contacts the patch. The unit will not make a noise or vibrate. Most people will not feel anything.




STEP 5 To increase the intensity turn the intensity wheel clockwise.



STEP 6 You can leave the patches on for up to 60 minutes allowing continued gently infusing. After removing patches, massage any residue into your skin.

INFUSION TREATMENT

We recommend Collagen and Retinol patches are only used once in 24 hours. If you have sensitive skin we would recommend leaving 48 hours between treatments.



STEP 1 Tear along the perforated edge of the Collagen and Retinol patch sachet. Take a patch from the plastic casing and peel off the liners to reveal the gel patch.



STEP 2 Place the patch onto the area to be treated. Place the other patch on the opposite side of the face.



STEP 3 Push the power button to ON. Set the micro-current button to (+)positive and set the intensity wheel to level 1



STEP 4 Slowly and gently move the roller back and forth on the patch.

Repeat the process over the other patch. Continue treatment by alternating between patches for 10-30 minutes (5-15 minutes on each patch).

NOTE: It is important that you hold the unit by the conductive panels. Correct operation is indicated by a blue light when the roller contacts the patch. The unit will not make a noise or vibrate. Most people will not feel anything.



STEP 5 To increase the intensity turn the intensity wheel clockwise.



STEP 6 You can leave the patches on for up to 60 minutes allowing continued gently infusing. After removing patches, massage any residue into your skin.